

## Xtalab Synergy-ED: To Cryo, Or Not to Cryo, That Is the ED Question

Dr. Mark Del Campo<sup>1</sup>, Dr. Khai-Nghi Truong<sup>2</sup>, Dr. Robert Bucker<sup>2</sup>, Dr. Fraser White<sup>2</sup>, Dr. Tomislav Stolar<sup>3</sup>, Dr. Krunoslav Užarević<sup>3</sup>, Dr. Luca Grisanti<sup>3</sup>, Dr. Ernest Meštrović<sup>4</sup>, Dr. Mathias Meyer<sup>5</sup>, Dr. Michał Jasnowski<sup>5</sup>, Dr. Akihito Yamano<sup>6</sup>, Dr. Sho Ito<sup>6</sup>, Dr. Eiji Okunishi<sup>7</sup>, Dr. Yoshitaka Aoyama<sup>7</sup>, Dr. Joseph Ferrara<sup>1</sup>  
<sup>1</sup>Rigaku Americas Corporation, <sup>2</sup>Rigaku Europe SE, <sup>3</sup>Ruder Bošković Institute, <sup>4</sup>University of Zagreb, <sup>5</sup>Rigaku Polska, <sup>6</sup>Rigaku Corporation, Haijima. <sup>7</sup>JEOL Ltd.  
[mark.delcampo@rigaku.com](mailto:mark.delcampo@rigaku.com)

The Rigaku XtaLAB Synergy-ED is a fully integrated electron diffractometer with a seamless workflow from data collection to 3D structure determination. The XtaLAB Synergy-ED is the result of Rigaku's collaboration with JEOL, synergistically combining each partner's core technologies: Rigaku's hybrid-pixel electron detector (HyPix-ED) and CrysAlisPro software, and JEOL's long-standing excellence in electron beam generation and control.

Using 3D electron diffraction (3DED), a.k.a. microED, single crystals of all classes below one- micron in size can be studied. The XtaLAB Synergy-ED offers a wide range of experiments, including low temperature, which confers benefits such as: reduced beam damage; improved resolution; reduced disorder; and in the case of electron diffraction, stabilization of samples in vacuo. As electron diffraction requires samples to be studied under high vacuum, cryo-transfer - freezing of samples prior to introduction to vacuum - is essential for sensitive compounds, e.g. proteins or MOFs. Combining cryo-transfer with variable detector distance further allows the study of solvated crystals, including proteins, in the XtaLAB Synergy-ED.

Here, we report two small molecule examples. Both are supramolecular solids featuring thymine (T) and 2,6-diaminopurine (D) in different stoichiometries (2:1 and 3:2). Thermal analyses (TGA-DSC) show that both samples are hydrated co-crystals. 3DED experiments at room temperature and high vacuum resulted in their anhydrous species (see Figure, left).

Cryo-transfer of flash-cooled crystals of the same batch measured at 100 K under high vacuum gave the 2:1 T-D monohydrate co-crystal and 3:2 T-D trihydrate co-crystal, respectively (see Figure, right). Calculated powder patterns from 3DED data match the experimental PXRD data, confirming phase purity. Our work unambiguously demonstrates that cryo-transfer techniques used with the XtaLAB Synergy-ED enabled the determination of hydrated structures, despite high vacuum conditions.



**Figure 1** Two examples illustrating cryo-transfer is important: (left) structures measured at ambient temperature and vacuum, and (right) using a cryo-transfer station at 100 K.

**Figure 1**