

MS3-O5 Automated solutions for fragment screening at the HZB MX beamlines

Karine M. Sparta¹, Uwe Mueller¹, Manfred S. Weiss¹, Monika Ühle¹, Franziska Huschmann¹, Johannes Schiebel², Andreas Heine², Gerhard Klebe²

1. Helmholtz-Zentrum Berlin, Macromolecular Crystallography, Berlin, Germany
2. University of Marburg, Department of Pharmaceutical Chemistry, Marburg, Germany

email: karine.sparta@helmholtz-berlin.de

X-ray diffraction plays an essential role in fragment-based lead discovery and design to identify target hits, and determine the location and structure of fragments bound to a target protein. Recently, a second beamline at the HZB MX facility [1] has been upgraded with a Pilatus detector and sample changer and is dedicated for fragment screening campaigns, enabling fast high throughput experiments to be performed routinely.

In order to help users cope with the huge amount of data thus created and reliably interpret their results rapidly, we have developed the expert software XDSAPP [2] for the analysis of diffraction images during measurements with minimal effort and time. It mainly uses the diffraction data processing program XDS [3], along with additional software like POINTLESS from the CCP4 suite [4], XDSSTAT [5], SFCHECK [6] and PHENIX.XTRIAGE [7] for automated decision making. An independent refinement pipeline for automated ligand search based on PHENIX and COOT [8] has been developed and can also be used within XDSAPP. In the presentation, the workflow of the pipeline will be described and its application to an example will be discussed.

XDSAPP is available free of charge for academic users from www.helmholtz-berlin.de/bessy-mx.

References

- [1] Mueller, U. et al. (2012). *J. Synchrotron Rad.* 19, 442-449.
- [2] Krug, M. et al. (2012). *J. Appl. Cryst.* 45, 568-572.
- [3] Kabsch, W. (2012). *Acta Cryst.* D66, 125-132.
- [4] Collaborative Computational Project (1994). *Acta Cryst.* D50, 760-763.
- [5] Diederichs, K. (2006). *Acta Cryst.* D62, 96-101.
- [6] Vaguine, A. et al. (1999). *Acta Cryst.* D55, 191-205.
- [7] Adams, P.D. et al. (2010). *Acta Cryst.* D66, 213-221.
- [8] Emsley, P. et al. (2010). *Acta Cryst.* D66, 486-501.

Keywords: XDSAPP, fragment screening, automation, software, data processing, refinement

MS4. Advances in phasing, refinement, and autobuilding

Chairs: Navraj Pannu, Eleonor Dodson

MS4-O1 ARCIMBOLDO, an *ab initio* approach to MR phasing

Isabel Uson¹, Massimo Sammito², Claudia Millán², Rafael J. Borges²

1. Institució Catalana de Recerca i Estudis Avançats (ICREA); Instituto de Biología Molecular de Barcelona (IBMB-CSIC)
2. Instituto de Biología Molecular de Barcelona (IBMB-CSIC)

email: uson@ibmb.csic.es

Ab Initio phasing of macromolecular structures, from the native intensities alone without need of derivatives or previous particular structural knowledge has been the object of a long quest. Currently, both experimental phasing and molecular replacement are going *ab initio*, as in the first case, inherent anomalous signal is succeeding with increasing generality[1] and in the second, the previous structural knowledge required is becoming unspecific[2,3,4]. Our own approach[5,6] relies on the combination of locating model fragments such as polyalanine alpha-helices, libraries of unspecific folds or portions of distant homologs with the program PHASER[7] and density modification with the program SHELXE[8]. Given the difficulties in discriminating correctly positioned fragments, many putative groups of fragments have to be tested in parallel, thus calculations are usually performed in a grid or supercomputer. The method has been called after the Italian painter Arcimboldo, who used to compose portraits out of fruits and vegetables. In the case of our program, most collections of fragments remain a "still-life", but some are correct enough for density modification to reveal the protein's true portrait. Latest developments, such as the single machine implementation ARCIMBOLDO_LITE, advances on the supercomputer GORDON at the SDSC or the first case of a solution for an unknown all-beta structure will be described. (http://chango.ibmb.csic.es/ARCIMBOLDO_LITE).

1) Liu, Guo, Chang, Cai, Assur, Mancina, Greene and Hendrickson. (2014). *Acta Cryst.* **D70**, 2544-2557.

2) DiMaio, Terwilliger, Read, Wlodawer, Oberdorfer, Wagner, Valkov, Alon, Fass, Axelrod, Das, Vorobiev, Iwai, Pokkuluri & Baker. (2011). *Nature*, **473**, 540-543.

3) Bibby, Keegan, Mayans, Winn & Rigden. (2012). *Acta Cryst.* **D68**, 1622-1631.

4) Shrestha & Zhang. (2015). *Acta Cryst.* **D71**, 304-312.

5) Millán, Sammito & Usón. (2015). *IUCrJ* **2**, 95-105.

6) Sammito, Millán, Rodríguez, de Iarduya, Meindl, De Marino, Petrillo, Buey, de Pereda, Zeth, Sheldrick & Usón. (2013). *Nature Methods*, **10**, 1099-1101

7) McCoy, Grosse-Kunstleve, Adams, Winn, Storoni, Read. (2007). *J. Appl. Cryst.* **40**, 658–674

8) Thorn & Sheldrick. (2013). *Acta Crystallogr. D* **69**, 2251-2256.

Keywords: Phasing, ab initio, macromolecular crystallography, fragment search, ARCIMBOLDO, BORGES, SHREDDER, PHASER, SHELXE

MS4-O2 ARP/wARP 7.5 and beyond

Victor S. Lamzin¹, Philipp Heuser¹, Joana Pereira¹, Nils Petersen¹

1. European Molecular Biology Laboratory, DESY,
Notkestrasse 85, 22607, Hamburg, Germany

email: victor@embl-hamburg.de

The ARP/wARP software for macromolecular model building is being continuously developed. Its recent release 7.5 (www.arp-warp.org), jointly with the CCP4 6.5, comes complete with additional innovations targeting protein model building at medium-to-low resolution range; including improved polypeptide recognition, NCS-restraints, atom update, estimation of solvent content and model accuracy, and an option to carry out SAD-refinement. The process of identifying and fitting of bound ligands now incorporates the 84 most common ligands and is able to use cif files defining bond, torsion and plane restraints.

Additionally, currently under development is the novel method of validating built protein fragments (presented in a dedicated poster at this meeting), as well as approaches to accurately estimating the quality of electron density maps and thus the phases. For the latter, a comprehensive scoring function has been designed that makes use of the distributions of third-order moment invariants computed throughout the density map. The distributions of the invariants have been shown to depend on the quality of phases, and the central moments of the distributions can conveniently capture the phase changes. The application of this method to decision-making, protocol guiding or validation in automated macromolecular structure solution and model building will be presented. This method is capable of ranking electron density maps within a resolution range from 2 to 6 Å and up to 80° phase error.

Keywords: ARP/wARP, software