

s5.m15.p2 **The crystal structure of a new mineral hillite $\text{Ca}_2(\text{Zn},\text{Mg})[\text{PO}_4]_2 \cdot 2\text{H}_2\text{O}$.** O. Yakubovich¹, W. Massa², R. Liferovich³, P. Gavrilenko¹, A. Bogdanova⁴, P. Tuisku⁵. ¹ Department of Crystallography, Moscow Lomonosov State University, Russia. ² Fachbereich Chemie der Philipps-Universität, Marburg, Germany. ³ Department of Geology, Lakehead University, Ontario, Canada. ⁴ Geological Institute, Kola Sc. Center, Apatity, Russia. ⁵ Institute of Geosciences, University of Oulu, Finland. E-mail: polinka_gav@mail.ru

Keywords: New mineral; Hillite; Crystal structure

Hillite, ideally $\text{Ca}_2(\text{Zn},\text{Mg})[\text{PO}_4]_2 \cdot 2\text{H}_2\text{O}$, is the zincian analogue of collinsite. Both hillite and collinsite belongs to the fairfieldite group, which includes phosphates and arsenates with general formula $A_2M[\text{TO}_4]_2 \cdot 2\text{H}_2\text{O}$, where $A = \text{Ca}, \text{Na}, \text{M} = \text{Mg}, \text{Fe}, \text{Zn}, \text{Ni}, \text{Co}, \text{Mn}$ and Cu , $T = \text{P}, \text{As}$. The new mineral was found in unmetamorphosed sediments of the Lower Cambrian Parachilna Formation at Reaphook Hill, South Australia. The single crystal of hillite do not exceed 50 μm in length; larger aggregates are inhomogeneous in terms of Zn:Mg ratio. The mineral is greenish and bluish to colourless, brittle with perfect cleavage along $\{010\}$ and $\{001\}$. It is biaxial positive and non-pleochroic; $D_{\text{meas}} = 3.16(2)\text{g}/\text{cm}^3$, $D_{\text{calc}} = 3.178\text{g}/\text{cm}^3$.

The crystal structure of hillite $\text{Ca}_2(\text{Zn}_{0.62}\text{Mg}_{0.38})[\text{PO}_4]_2 \cdot 2\text{H}_2\text{O}$, triclinic, space group PI , a 5.736(1), b 6.767(2), c 5.462(1)Å, α 97.41(2), β 108.59(2), γ 107.19(2)°, $V = 186.05(8)\text{Å}^3$, $Z = 1$, has been determined [automated single-crystal diffractometer, $\text{MoK}\alpha$, graphite monochromator, imaging-plate area detector system, $T = 293\text{ K}$, 2928 reflections, $wR_2 = 0.0998$ for all 1078 unique reflections, $R = 0.0378$ for 993 observed reflections with $I \geq 2\sigma(I)$]. The refinement of site occupancies showed that Zn and Mg are in the octahedral position at the centre of symmetry in the proportion 0.623(5):0.377(5). The positions of two independent H atoms were obtained from difference-Fourier syntheses and were refined under isotropic approximation. The geometric characteristics of hydrogen bonds are determined, and bond-valence analysis is made. Interatomic distances in hillite structure were analysed and compared with those of other members of the fairfieldite group.

The crystal structure of hillite consists of isolated $\text{MO}_4(\text{H}_2\text{O})_2$ octahedra, which are connected by orthophosphate tetrahedra to form chains parallel to the c axis. Each tetrahedron shares two vertices with neighbouring octahedral along the chain, whereas the two other oxygen vertices co-ordinate a Ca^{2+} ion. In the a and b directions, these chains are held together by Ca atoms and hydrogen bonds.

According to the geometry of the hydrogen bonds the fairfieldite group of minerals can be divided into two subgroups (fairfieldite and collinsite). We consider the new mineral hillite as a member of the collinsite subgroup.

s5.m15.p3 **Birefringence imaging of minerals using a tilting stage.** A.M. Glazer & L. Pajdzik, *Physics Department, University of Oxford, Parks Rd., Oxford OX1 3PU, UK. E-mail: glazer@physics.ox.ac.uk*

Keywords: Birefringence; Minerals; Microscopy

We show that by combining the Metripol imaging system (see www.metripol.com) with a tilting stage it is possible to obtain very precise birefringence information on crystals. The Metripol microscope uses a combination of a rotating polarizer and a circular analyser to separate out three types of image, one representing the light transmission through the specimen, one showing the orientation of the optical indicatrix at any point in the image, and one giving quantitative information on $|\text{sind}|$, where d is the phase shift of the light. The data are collected as a function of two angles of tilt of the microscope stage it is possible to solve equations to give birefringence information both in-plane and out of plane of the specimen. This method shows potential for the use of automatically identifying crystals, especially minerals in microscope rock sections.