

MS10-O2 Molecular Interactions in Pharmaceutical CompoundsKatharina Edkins¹

1. School of Medicine, Pharmacy and Health, Durham University, United Kingdom

email: katharina.fucke@durham.ac.uk

Small organic molecules, especially in the pharmaceutical sciences, tend to crystallise in a plethora of different crystal forms, either as pure compounds or with the inclusion of solvent molecules. Due to their changed physico-chemical characteristics, such as melting point, compressibility, solubility and thus bioavailability, and physical and chemical stability, different crystal forms can pose a problem to the manufacture of medicines.^[1] It is thus crucial to understand the crystallisation behaviour and manufacturability of these compounds in order to avoid problems in the life-time of the medicine and costly recalls comparable to ritonavir^[2] or rotigotine.^[3] Bioactive molecules and pharmaceuticals typically have multiple functional groups, enabling them to interact with receptors and thus show pharmacological action. In the solid-state, the interactions through these functional groups are the driving forces of molecular recognition. By applying X-ray and neutron diffraction methods as well as thermoanalysis, vapour sorption and spectroscopic analysis in combination with computational techniques, we are probing the strong and weak interactions within the crystal forms and during the crystallisation in order to understand and predict their characteristics.

[1] R. Hilfiker, *Polymorphism: In the Pharmaceutical Industry*, Wiley-VCH Verlag GmbH & Co. KGaA, Weinheim, Germany, **2006**. [2] S. R. Chemburkar, J. Bauer, K. Deming, H. Spiwek, K. Patel, J. Morris, R. Henry, S. Spanton, W. Dziki, W. Porter, J. Quick, P. Bauer, J. Donaubaue, B. A. Narayanan, M. Soldani, D. Riley and K. McFarland, *Organic Process Research & Development* **2000**, *4*, 413-417. [3] I. B. Rietveld and R. Ceolin, *Journal of Pharmaceutical Sciences* **2015**, *104*, 4117-4122.

Keywords: Polymorphism, Pharmaceutical, Molecular recognition**MS10-O3** High resolution neutron and X-ray diffraction RT studies of an H-FABP – Oleic acid complex: study of the internal water cluster and the ligand binding by a transferred multipolar electron density distributionAlberto Podjarny¹, Eduardo I. Howard², Benoit Guillot³, Matthew P. Blakeley⁴, Michael Haertlein⁵, Martine Moulin⁵, Andre Mitscher⁶, Alexandra Cousido-Siah¹, Firas Fadel¹, Wanda Valsecchi⁹, Takashi Tomizaki⁷, Tania Petrova⁸, Julien Claudot³

1. IGBMC, CNRS, INSERM, UdS, Illkirch, France
2. IFLYSIB, CONICET, UNLP, La Plata, Argentina
3. CRM2, CNRS, U.Lorraine, Vandoeuvre-les-Nancy, France
4. ILL, Grenoble, France
5. PBS, Grenoble, France
6. IQFB, UBA, Buenos Aires, Argentina
7. SLS, Villigen PSI, Switzerland
8. IMPB, RAS, Pouchino, Russia

email: apodjarny@gmail.com

Crystal diffraction data of heart fatty acid binding protein (H-FABP) in complex with oleic acid were measured at room temperature with high resolution X-ray and neutron protein crystallography (0.98 Å and 1.90 Å resolution respectively). These data provided very detailed information about the cluster of water molecules and the bound oleic acid in the H-FABP large internal cavity. The jointly refined X-ray/neutron structure of H-FABP was complemented by a transferred multipolar electron density distribution using the parameters of the ELMAM2 library. The resulting electron density allowed a precise determination of the electrostatic potential in the fatty acid (FA) binding pocket. Bader's quantum theory of atoms in molecules was then used to study interactions involving the internal water molecules, the FA and the protein. This approach showed H-H contacts of the FA with highly conserved hydrophobic residues known to play a role in the stabilization of long chain FAs in the binding cavity. The determination of water hydrogen (deuterium) positions allowed the analysis of the orientation and electrostatic properties of the water molecules in the very ordered cluster. As a result, a significant alignment of the water molecules permanent dipoles with the protein electrostatic field was observed. This can be related to the dielectric properties of hydration layers around proteins, where the shielding of electrostatic interactions depends directly on the rotational degrees of freedom of the water molecules in the interface.

Keywords: Neutron protein crystallography; High resolution room temperature X-ray crystallography; Fatty acid binding protein; Protein hydration layer; AIM Topological properties